

# Inferring baseline optical properties of the human head

Alex Barnett

`barnett@nmr.mgh.harvard.edu`

MGH/MIT/HMS AAM NMR CBI

# Aims

- How well can we measure baseline optical tissue properties? Given...
  - 3d anatomical MRI data
  - optically-uniform segmented tissue types
  - time-resolved measurements
  - single optical  $\lambda$
- Motivations:
  - functional imaging requires accurate baseline properties
  - more  $\lambda$ 's  $\rightarrow$  absolute [Hb] and [HbO].
  - sets an upper bound on capability *without* MRI data.

# Outline

1. Bayesian method overview
2. simple layer system
3. likelihood
4. results in layer
5. optode calibration & location
6. preliminary head
7. issues & conclusion

# Method overview


What do measurements  $y$  tell you about parameters  $x$  ?

Inference  $\rightarrow$  probability distribution functions (PDFs)

# Method overview

What do measurements  $y$  tell you about parameters  $x$  ?

Inference  $\rightarrow$  probability distribution functions (PDFs)




$$\begin{aligned} p(\mathbf{x}|\mathbf{y}) &\propto p(\mathbf{y}, \mathbf{x}) \\ &= p(\mathbf{y}|\mathbf{x}) \cdot p(\mathbf{x}) \\ \text{posterior} &\propto \text{likelihood} \cdot \text{prior} \end{aligned}$$

Constant prior  $\Rightarrow$  look for peaks in likelihood.

# Method overview

What do measurements  $y$  tell you about parameters  $x$  ?

Inference  $\rightarrow$  probability distribution functions (PDFs)



$$\begin{aligned} p(\mathbf{x}|\mathbf{y}) &\propto p(\mathbf{y}, \mathbf{x}) \\ &= p(\mathbf{y}|\mathbf{x}) \cdot p(\mathbf{x}) \\ \text{posterior} &\propto \text{likelihood} \cdot \text{prior} \end{aligned}$$


Constant prior  $\Rightarrow$  look for peaks in likelihood.

- No artificial regularization
- Peak widths give all errorbars & correlations

# Method overview

What do measurements  $y$  tell you about parameters  $x$  ?

Inference  $\rightarrow$  probability distribution functions (PDFs)




$$\begin{aligned} p(\mathbf{x}|\mathbf{y}) &\propto p(\mathbf{y}, \mathbf{x}) \\ &= p(\mathbf{y}|\mathbf{x}) \cdot p(\mathbf{x}) \\ \text{posterior} &\propto \text{likelihood} \cdot \text{prior} \end{aligned}$$

Constant prior  $\Rightarrow$  look for peaks in likelihood.

- No artificial regularization
- Peak widths give all errorbars & correlations

Currently: testing with numerically-generated noisy measurements  $y$

# Simple 2-layer system




S-D separations of 7, 14, 21, 28 mm

Parameter vector  $\mathbf{x} \equiv [\mu_a(1), \mu_a(2), \mu'_s(1), \mu'_s(2)]$



# Likelihood

- $f(\mathbf{x}) = \text{forward model}$  (signal expectation)
- $p_{\text{noise}} = \text{noise model}$




$$p(\mathbf{y}|\mathbf{x}) = p_{\text{noise}}(\mathbf{y} | \mathbf{f}(\mathbf{x}))$$


uncorr. gaussian  $\longrightarrow \prod_m \frac{1}{\sqrt{2\pi}\sigma_m(\mathbf{x})} e^{-\frac{1}{2} \frac{[y_m - f_m(\mathbf{x})]^2}{\sigma_m^2(\mathbf{x})}}$

$\sigma$  is some (growing) function of  $\mathbf{f}$ , giving detection statistics.

# Look at sensitivity




# Sensitivity compared to noise



$\sigma$ -normalized changes :  $\beta_m \equiv \Delta f_m / \sigma_m$

# Maximizing Likelihood




Minimize 'objective function'  $\text{NLL} \equiv -\ln p(y|x)$

- gaussian noise  $\rightarrow \approx$  'weighted least squares'
- peak very narrow in  $x_{\text{layer 1}}$   $\rightarrow$  I show only  $x_{\text{layer 2}}$
- 1-2 minutes per optimization


# Results : photon number

typ tissue properties  $\mathbf{x}_{\text{expt}} = (0.01, 0.01, 2, 1) \text{ mm}^{-1}$




- more photons  $\rightarrow$  narrower peak
- true  $\mathbf{x}_{\text{expt}}$  rarely outside peak — good!

# Results : varying $\mu_a(2)$



- other 3 parameters held constant
- photon # : 67 photons at det4
- realistic inference of errorbars

# Results : varying $\mu'_s(2)$



Generally good agreement. Reliability problems...

- noise model mismatch ? / optimization getting stuck


# Integrating out free parameters

Width in  $\mathbf{x}_{\text{layer 1}}$  is *much* less than in  $\mathbf{x}_{\text{layer 2}}$ .

We only care about  $\mathbf{x}_{\text{layer 2}}$  (e.g. cortex in head).

Once peak found, use gaussian approx: analytic integral over

$\mathbf{x}_{\text{layer 1}}$ :



$$\int d\mathbf{x} e^{-\frac{1}{2}\mathbf{x}^T H \mathbf{x}} = \frac{(2\pi)^{N/2}}{(\det H)^{1/2}}$$

This illustrates the general Bayesian recipe for free parameters :  
integrate over them.



# Optode calibration & placement

*Optode calibration* :  $(N_s + N_d - 1)$  free scale parameters

- As for layer 1, they will be narrow-width
- integrate out with gaussians (fast)


*Placement* : choose best source/detector locations

- use peak volume  $(\det H)^{-1/2}$  as objective func.
- fix  $\mathbf{x} = \mathbf{x}_{\text{expt}}$ , and optimize over locations.

For gaussian noise model  $H \approx J^T \cdot \text{diag}(1/\sigma) \cdot J$ .  
with jacobian  $J_{mn}(\mathbf{x}) \equiv \partial f_m / \partial x_n$ .

# Noise model details

Used uncorrelated gaussian model:




- gaussian approx to poisson, clipped at both ends
- Collect more photons  $\Rightarrow$  model error dominates
- Other more *robust* noise models (power law tails, etc) possible, easy to implement in Bayesian formalism.

# Forward model details

Time-resolved detector signals  $f$  given params  $x$ .

Written finite-difference time-domain (FDTD) code:



- arbitrary 3d tissue geometries
- 0.5s per source, small system  $6\text{cm} \times 6\text{cm} \times 3\text{cm} \times 2\text{ns}$
- Diffusion Approx, validated against Monte Carlo
- Robin BCs, surface normals only  $\pm xyz$ .
- evolution: 'forward-Euler'  $O(\Delta t)$ , small  $\mu'_s$  slows it down.


# Forward model issues

There are  $O(\Delta t^2)$  methods ('implicit', e.g. ADI) :

- faster (less timesteps), but nonsmooth fluence *bad!*

## Boundary Conditions

- *do* matter.
- 'Stiffness' tricky for FDTD stability



Avoid large system (head) by matching to  $\infty$  :

fluence components  $\omega \ll c\mu_a$  obey Helmholtz eqn with *fixed*  $k \approx i\sqrt{3\mu_a\mu'_s}$ . So, 'radiative' BC is just Robin BC.


# Nonlinear optimization issues

Forward model with discontinuities (jumps) = *bad* :

ridges in  $f(\mathbf{x})$

fake local minima

Had to be removed!



Derivative info *vastly* improves speed/robustness:

- Adjoint ('reverse') differentiation: get  $\nabla_{\mathbf{x}} f$  wrt *all*  $x_n$  with little more effort than  $f$  (e.g. Hielscher, Klose, Hanson 1999)

# Future directions

- MRI-segmented heads: 5 tissue-types, effect of CSF.

# Future directions

- MRI-segmented heads: 5 tissue-types, effect of CSF.
- Forward model improvements: adjoint differentiation, ADI timestepping.

# Future directions

- MRI-segmented heads: 5 tissue-types, effect of CSF.
- Forward model improvements: adjoint differentiation, ADI timestepping.
- Use peak width (Hessian) to optimize optode locations



# Future directions

- MRI-segmented heads: 5 tissue-types, effect of CSF.
- Forward model improvements: adjoint differentiation, ADI timestepping.
- Use peak width (Hessian) to optimize optode locations
- Test integrating out free params: optode calibration, start time...

# Future directions

- MRI-segmented heads: 5 tissue-types, effect of CSF.
- Forward model improvements: adjoint differentiation, ADI timestepping.
- Use peak width (Hessian) to optimize optode locations
- Test integrating out free params: optode calibration, start time...
- Noise models: Poisson for low photon counts, detector saturation.

# Future directions

- MRI-segmented heads: 5 tissue-types, effect of CSF.
- Forward model improvements: adjoint differentiation, ADI timestepping.
- Use peak width (Hessian) to optimize optode locations
- Test integrating out free params: optode calibration, start time...
- Noise models: Poisson for low photon counts, detector saturation.
- More sources, experimental phantom verification, heads...

# Conclusions

- Time-resolved data with few detected photons can infer optical parameters in layer 8mm below surface, to  $\pm 10\%$

# Conclusions

- Time-resolved data with few detected photons can infer optical parameters in layer 8mm below surface, to  $\pm 10\%$
- Full posterior (errorbars & correlations) can be handled

# Conclusions

- Time-resolved data with few detected photons can infer optical parameters in layer 8mm below surface, to  $\pm 10\%$
- Full posterior (errorbars & correlations) can be handled
- Developed & validated rapid 3d diffusion forward model

# Conclusions

- Time-resolved data with few detected photons can infer optical parameters in layer 8mm below surface, to  $\pm 10\%$
- Full posterior (errorbars & correlations) can be handled
- Developed & validated rapid 3d diffusion forward model
- Bayesian optode calibration and optimal location recipes